


MEDICAL LIGHT GUIDES


Input Fittings

1 LETTER




Distal Fittings

2 LETTERS



Specifying Information

6 position part number describes



1. Three Choices: _____
 - *MLG - Standard Medical Light guide
 - *MMLG - Micro Medical Light Guide
 - **UMLG - Universal Medical Light guide
2. ILLUMINATOR END FITTING (1 letter) _____
3. BUNDLE DIAMETER (2mm-6mm) _____
4. LENGTH (6, 7.5 or 10Ft) _____
5. DISTAL END (2 Letters) _____
6. COLOR - Pick from the list _____

1 2 3 4 5 6
 MLG S 5.0 7.5 ST GRAY

- COLOR**
- Clear
 - Red
 - Orange
 - Yellow
 - Pink
 - Blue
 - Purple
 - Green
 - Black
 - Light blue
 - Gray

*When specifying MLG - select one Input and Distal end.
 **When specifying UMLG - Any combination may be used.